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For HEMOCHRON® Jr. Microcoagulation Systems

Package Insert
English

INTENDED USE
The directCHECK® Whole Blood Quality Controls are dried whole blood preparations which have been assayed and are intended to be used to perform quality control assays using the HEMOCHRON Jr. Microcoagulation test cuvettes.

For in vitro Diagnostic Use

SUMMARY AND EXPLANATION
Blood coagulation instruments and assays should be quality controlled prior to and during routine use. Performance ranges are provided with each control product against which the user should compare results. A proper quality assurance program includes regular instrument maintenance, quality control assay measurements and complete performance records.

REAGENTS
Level 1 and Level 2 quality control products are provided separately. These preparations consist of dried fixed bovine red blood cells, buffered sheep and horse plasma. Assayed clotting time values are provided with each lot of material.

Each control preparation is provided in a dropper vial. Each dropper vial also contains diluent used to rehydrate the dried whole blood control. These diluent preparations consist of distilled water, sodium chloride, Tween® 20, ProClim®, anticoagulant and calcium chloride. Calcium chloride is used for recalcification of the citrated whole blood preparation.

Four Level 1 quality control products are provided separately, to be used with the following HEMOCHRON Jr. tests:

- APTT Level 1 quality control (cat. # DCJAPTT-N); for use with HEMOCHRON Jr. APTT test cuvette in any model HEMOCHRON Jr. instrument except original HEMOCHRON Jr.
- PT Level 1 quality control (cat. # DCJPT-N) for use with HEMOCHRON Jr. PT test cuvette in any model HEMOCHRON Jr. instrument.
- ACT+ Level 1 quality control (cat. # DCJACT-A) for use with HEMOCHRON Jr. ACT+ test cuvette in any model HEMOCHRON Jr. instrument.
- ACT-LR Level 1 quality control (cat. #DCJLR-N) for use with HEMOCHRON Jr. ACT-LR test cuvette in any model HEMOCHRON Jr. instrument.

Four Level 2 quality control products are provided separately, to be used with the following HEMOCHRON Jr. tests:

- APTT Level 2 quality control (cat. # DCJAPTT-A); for use with HEMOCHRON Jr. APTT test cuvette in any model HEMOCHRON Jr. instrument except original HEMOCHRON Jr.
- PT Level 2 quality control (cat. # DCJPT-A); for use with HEMOCHRON Jr. PT test cuvette in any model HEMOCHRON Jr. instrument.
- ACT+ Level 2 quality control (cat. # DCJACT-A); for use with HEMOCHRON Jr. ACT+ test cuvettes in any model HEMOCHRON Jr. instrument.
- ACT-LR Level 2 quality control (cat. # DCJLR-A); for use with HEMOCHRON Jr. ACT-LR test cuvettes in any model HEMOCHRON Jr. instrument.

WARNINGS AND PRECAUTIONS
CAUTION: Animal blood products are included in the directCHECK quality control preparation. No material in any preparation is derived from human blood. All blood products, whether of human origin or animal origin, however, should be handled with care.

CAUTION: All used test cuvettes and directCHECK vials should be considered as potentially infectious, handled with care and disposed of in accordance with standard medical waste disposal policy.

HANDLING AND STORAGE
When refrigerated (2-8°C) the vials are stable until the marked expiration date. The quality control product should never be exposed to temperatures in excess of 37°C. Reconstituted vials should be used immediately. directCHECK Quality Control products may also be stored at room temperature for up to 4 weeks. (The marked expiration date must not be exceeded.) A retraining label is provided and should be marked with 4 weeks dating if room temperature storage is selected.

PROCEDURE
Materials Provided
directCHECK quality control product kits for use with the HEMOCHRON Jr. Microcoagulation Systems contain the following:

- 15 dropper vials of dried whole blood control (0.5 mL) in a glass ampule provided with 0.7 mL of diluent.
- 4 reusable protective sleeves, for use in crushing ampules.

TROUBLESHOOTING
Please refer to the appropriate instrument Operator's Manual "Troubleshooting" section for instrument software fault messages, definitions and actions.

PROBLEM	CAUSE	ACTION
Control value is below published range	Reconstituted control was not thoroughly mixed.	Repeat test, ensuring that inner glass ampule is crushed at least two times and reconstituted control is thoroughly mixed by inverting the vial end to end 10 times.
	Time period between control material mixing and addition to test cuvette is too long.	Repeat test. After proper mixing of control material and diluent, IMMEDIATELY discard first drop and then dispense reconstituted control material into cuvette well. Do not remove protective sleeve until after control material has been dispensed from vial into cuvette well. Ensure vial label is removed prior to crushing vial.
Control value is above published range	Inner glass ampule was not adequately crushed.	Repeat test, ensuring that inner glass ampule is crushed at least two times prior to inverting the vial end to end 10 times.
	Reconstituted control was not thoroughly mixed.	Repeat test, ensuring that reconstituted control is thoroughly mixed by inverting the vial end to end 10 times.
	Vial cap is removed prior to inverting, allowing diluent to leak from vial.	Repeat test, ensuring that vial cap is not removed prior to inverting. Remove vial cap when control material is to be dispensed into cuvette well.
Fault message obtained (see appropriate HEMOCHRON Jr. Operator's Manual)	Reconstituted control was not thoroughly mixed.	Repeat test, ensuring that reconstituted control is thoroughly mixed by inverting the vial end to end 10 times.
	Reconstituted control was not thoroughly mixed.	Repeat test, ensuring that inner glass ampule is crushed at least two times prior to inverting end to end 10 times.
	Bubbles are present in control material when added to cuvette well.	Repeat test, avoiding excessively vigorous shaking that may cause bubbles or foaming.
		Repeat test, ensuring reconstituted control is thoroughly mixed by inverting the vial end to end 10 times.

IF OUT-OF-RANGE RESULTS OR FAULT MESSAGES PERSIST, CONTACT ACCRIVA TECHNICAL SUPPORT.

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