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For HEMOCHRON® Jr.
Microcoagulation Systems

Package Insert
English

INTENDED USE

The directCHECK® Citrate Whole Blood Quality Controls are dried whole blood preparations which have been assayed and are intended to be used to perform quality control assays using the HEMOCHRON® Jr. Microcoagulation Citrate PT and Citrate APTT test cuvettes.

For in vitro Diagnostic Use

SUMMARY AND EXPLANATION

Quality control instruments and assays should be quality controlled prior to and during routine use. Performance ranges are provided with each Control Product against which the user should compare results. A proper quality assurance program includes regular instrument maintenance, quality control assay measurements and complete performance records.

REAGENTS

Level 1 and Level 2 quality control products are provided separately. These preparations consist of dried fixed bovine red blood cells, buffered sheep and horse plasma. Assayed clotting time values are provided with each lot of material.

Each control preparation is provided in a dropper vial. Each dropper vial also contains diluent water, sodium chloride, Tween 20, and anticoagulant.

Level 1 Citrate PT and Citrate APTT quality control (DQCPFN) is to be used with the HEMOCHRON Jr. Citrate PT and Citrate APTT test using a HEMOCHRON Jr. II or Jr. Signature instrument.

The Level 2 Citrate PT quality control (DQCPFA) is provided separately, to be used with the Citrate PT test using a HEMOCHRON Jr. II or Jr. Signature instrument.

The Level 2 Citrate APTT quality control (DQCAPTA) is provided separately, to be used with the Citrate APTT test using a HEMOCHRON Jr. II or Jr. Signature instrument.

CAUTION: Animal blood products are included in the directCHECK quality control preparation. No material in any preparation is derived from human blood. All blood products, whether of human origin or animal origin, however, should be handled with care.

HANDLING AND STORAGE

When refrigerated (2–8°C) the vials are stable until the marked expiration date. The quality control product should never be exposed to temperatures in excess of 37°C. Reconstituted vials should be used immediately. directCHECK Quality Control products may also be stored at room temperature for up to 4 weeks. (The marked expiration date must not be exceeded.) A retarding label is provided and should be marked with 4 weeks dating if room temperature is selected.

PROCEDURE

Materials Provided

directCHECK Citrate quality control product kits for use with the HEMOCHRON Jr. II and Jr. Signature Microcoagulation Systems contain the following:

- 15 dropper vials of dried whole blood control (0.5 ml) in a glass ampule provided with 0.7 ml of diluent.
- 4 reusable protective sleeves, for use in crushing ampules.

Materials Required But Not Provided

- HEMOCHRON Jr. or Signature instrument
- HEMOCHRON Jr. Citrate PT test cuvettes
- HEMOCHRON Jr. Citrate APTT test cuvettes

Preparation of Control Material

A. HEMOCHRON Jr. Assays
Remove the appropriate test cuvettes from the refrigerator and allow them to come to room temperature (15–30°C) prior to testing. This should require up to 60 minutes.

B. Quality Control Product

Remove the appropriate quality control dropper vials from the refrigerator and allow them to come to room temperature. This should require up to 60 minutes. Visually inspect vial to ensure that the glass ampule is intact.

Test Procedure

CAUTION: All used test cuvettes and directCHECK vials should be considered as potentially infectious, handled with care and disposed of in accordance with standard medical waste disposal policy.

- Insert cuvette (room temperature) into the cuvette slot on the side of instrument.
- During pre-warm stage, observe display for fault/warning messages.
- The instrument will signal ready with an audible tone, the display will indicate "Add Sample" and "Press Start".
- Reconstitute the (room temperature) dropper vial contents as follows:

NOTE: Reconstitution and mixing of the whole blood control material should be accomplished quickly and without delay in any step. Once the dried control material has been reconstituted, the sample should be used immediately, as clotting will occur.

Remove the label from the vial. Insert vial into protective sleeve. Holding vial upright, tap the vial on a table top to settle glass ampule to the bottom of the vial. Crush the inner glass ampule by either bending the vial over the edge of a table top or by crushing the vial between two fingers. Immediately repeat the crushing action one to two additional times to ensure complete breakage of the glass ampule.

Quickly invert the dropper vial end to end 10 times.

5. While inverting the vial (dropper tip down), use a downward snapping motion of the wrist to ensure the control material flows to dropper tip. Remove and retain the vial cap. Squeeze the vial to discard the first drop of the control material. Immediately dispense as many drops of control material as needed to fill the cuvette sample well flush to the top. Should a large dome extend over the top of the center sample well, push it over into the outer sample well.

6. Depress the START key.

7. Recap the control vial and remove the control vial from the protective sleeve. Dispose of the vial and vial cap appropriately and retain the protective sleeve for reuse.

8. Wait for the single beep signaling end of the test. Two beeps indicate a fault condition.

9. Record results.

NOTE: Refer to the "Troubleshooting" section should any fault message appear on the display during this procedure.

PERFORMANCE CHARACTERISTICS

Expected Values

Acceptable ranges for each Quality Control Product are published on the back of this package insert. These values have been established at the manufacturer's facility and are lot specific. The range represents the acceptable performance ranges for all HEMOCHRON Jr. II and Jr. Signature instruments. The user should be able to obtain a value within the published range included with each product. A proper quality assurance program includes regular instrument maintenance, quality control measurements and maintenance of complete performance records.

ITC recommends that each institution establish its own expected range of response based on the mean ±2 standard deviations of at least 20 repeated test results. The local mean values established should fall within the manufacturer's acceptable performance range. Studies show that intra-laboratory variation in test results should produce a coefficient of variation of approximately 14% or less for coagulation control tests.

LIMITATIONS

Technique variation can significantly alter performance characteristics. While the HEMOCHRON Jr. Systems have a limited number of operator variables that can affect performance, improper handling of the Quality Control Product can significantly affect results. Care must be taken to establish and adhere to the recommended test procedure found in this package insert and the associated test cuvette package insert. Variations in ambient temperature should be minimized for consistency of test performance. Avoid vigorous shaking which may cause bubbles or foaming. In cases where results fall outside of the assigned range, patient results should be considered suspect and performance of quality control tests should be repeated.

Published performance ranges reflect results found in most clinical settings. Each institution should establish its own expected range of results and provide a data tracking record to identify instruments, cuvettes, reagents or operators whose performance is outside of the institution's expected range.

In cases where results are outside of the expected/published range, the cause will likely fall into one of the following possible categories.

- Operator technique
- Quality Control Product / test cuvette
- Instrumentation
- Environmental Factors

TROUBLESHOOTING

Please refer to the appropriate instrument Operator's Manual "Troubleshooting" section for instrument software fault messages, definitions and actions.

PROBLEM	CAUSE	ACTION
Control value is below published range	Reconstituted control was not thoroughly mixed.	Repeat test, ensuring that inner glass ampule is crushed at least two times and reconstituted control is thoroughly mixed by inverting the vial end to end 10 times.
	Time period between control material mixing and addition to test cuvette is too long.	Repeat test. After proper mixing of control material and diluent, IMMEDIATELY discard first drop and then dispense reconstituted control material into cuvette well. Do not remove protective sleeve until after control material has been dispensed from vial into cuvette well. Ensure vial label is removed prior to crushing vial.
Control value is above published range	Inner glass ampule was not adequately crushed.	Repeat test, ensuring that inner glass ampule is crushed at least two times prior to inverting the vial end to end 10 times.
	Reconstituted control was not thoroughly mixed.	Repeat test, ensuring that reconstituted control is thoroughly mixed by inverting the vial end to end 10 times.
	Vial cap is removed prior to inverting, allowing diluent to leak from vial.	Repeat test, ensuring that vial cap is not removed prior to inverting. Remove vial cap when control material is to be dispensed into cuvette well.
Fault message obtained (see appropriate HEMOCHRON Jr. Operator's Manual)	Reconstituted control was not thoroughly mixed.	Repeat test, ensuring that reconstituted control is thoroughly mixed by inverting the vial end to end 10 times.
	Bubbles are present in control material when added to cuvette well.	Repeat test, avoiding excessively vigorous shaking that may cause bubbles or foaming.
		Repeat test, ensuring reconstituted control is thoroughly mixed by inverting the vial end to end 10 times.
		Repeat test, ensuring that reconstituted control material flows to dropper tip by using a downward snapping motion of the wrist prior to dispensing into cuvette sample well.

IF OUT-OF-RANGE RESULTS OR FAULT MESSAGES PERSIST, CONTACT ITC TECHNICAL SERVICES.

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directCHECK® and ProTime® are registered trademarks of ITC.