# HEMOCHRON®

## Whole Blood Coagulation Systems

### Activated Clotting Time (ACT) FTCA510/HRFTCA510, FTK-ACT/HRFTK-ACT & P214

## Package Insert

English

#### INTENDED USE

The HEMOCHRON® Activated Clotting Time (ACT) is commonly used for heparin anticoagulation monitoring during cardiopulmonary bypass surgery, percutaneous transluminal coronary angioplasty (PTCA), interventional radiology, extracorporeal membrane oxygenation (ECMO), hemofiltration, hemodialysis and critical care.

The various ACT tests are intended solely for use with HEMOCHRON models Response, 8000, 801 and 401. The tests are performed using fresh whole blood at the patient bedside.

For *in vitro* Diagnostic Use

#### SUMMARY AND EXPLANATION

Beside coagulation systems eliminate many variables known to affect the reliability of activated clotting time testing, including specimen collection, handling and processing. HEMOCHRON ACT tests are sensitive to different levels of heparin anticoagulation, making them useful in a variety of clinical settings.

Whenever blood is exposed to a foreign surface such as an extracorporeal circuit, the clotting process is triggered. In especially invasive procedures like bypass surgery and PTCA, moderate to high levels of heparin anticoagulation are required to prevent thrombosis. In procedures such as ECMO and hemodialysis, more moderate levels of heparin anticoagulation are necessary. Low level heparin anticoagulation may also be administered concurrently with thrombolytic therapy, prophylactically after myocardial infarction and some surgeries to prevent the progression of thrombosis, and in the treatment of deep vein thrombosis, and pulmonary embolism.

The ACT test, first described by Hattershy in 1966,<sup>1</sup> is the method of choice for monitoring heparin therapy pre-, intra-, and post-bypass and during PTCA.<sup>2</sup>

While heparin therapy is essential in maintaining hemostasis during these procedures, its administration poses significant risk to the patient. Heparin is said to vary twelve-fold in both patient sensitivity and the rate at which it disappears from the blood!<sup>3</sup> Overdosing heparin can result in dangerous bleeding, whereas underdosing heparin can lead to thrombosis. Therefore, monitoring heparin therapy is vital in guarding against these undesirable side effects.

The ACT is a test in which fresh whole blood is added to a test tube containing an activator (Celite®, glass particles or kaolin) and timed for the formation of a clot. The type of activator used may affect the degree of prolongation of the ACT.

Celite® (diatomaceous earth) has long been hailed as the standard for high level heparin monitoring due to its excellent activating properties. However, protease inhibitors such as aprotinin, which may be administered to reduce post-operative bleeding, especially during cardiopulmonary bypass surgery, can prolong the Celite-activated ACT. Kaolin is unaffected by moderate doses of aprotinin. The kaolin-activated ACT can be used as an alternative ACT for heparin monitoring in patients receiving this drug.<sup>4</sup> ACT tubes containing activators like glass particles are most appropriate when low levels of heparin are used.<sup>4</sup>

The HEMOCHRON coagulation system offers a full range of ACT tubes for every level of heparin anticoagulation, from prophylaxis to intensive. The test may be performed at bedside where the blood sample is drawn directly into a test tube or is transferred into one by syringe. That test tube is then inserted into the instrument. Mechanical detection of a fibrin clot in the blood sample automatically terminates the test, causing the instrument's digital timer to display the coagulation time in seconds.

#### REAGENTS

##### ACT for high heparin level monitoring

*Bypass, Vascular Surgery and PTCA*

(FTCA510/HRFTCA510 - bar-coded for Response)

- Flip-top non-evacuated glass test tube with diatomaceous earth

*Heparin surgery in the presence of Aprotinin*

(FTK-ACT/HRFTK-ACT - bar-coded for Response)

- Flip-top non-evacuated glass test tube with kaolin

##### ACT for low to moderate heparin level monitoring

*Hemodialysis, Hemofiltration and ECMO*

(P214)

- Flip-top plastic test tube with glass beads

*NOTE: Information stated in this package insert for the Celite and kaolin ACT tubes for test procedure, performance characteristics, limitations, etc., are identical for both the bar-code labeled tubes (HRFTCA510 and HRFTK-ACT) and the non bar-code labeled tubes (FTCA510 and FTK-ACT).*

#### STORAGE AND STABILITY

ACT tests are stable when stored at room temperature and used prior to the expiration date on the carton and case.

#### SPECIMEN COLLECTION AND HANDLING

##### Materials Provided

- Activated Clotting Time test tubes
  - (FTCA510/HRFTCA510) black flip-top
  - (FTK-ACT/HRFTK-ACT) gold flip-top
  - (P214) clear flip-top

- HEMOCHRON Response, 8000, 801 or 401
- Appropriate syringe for sample collection

Use a 5cc syringe. Adhere to the appropriate technique (A, B or C).

#### Indwelling venous blood-line

- Do not obtain blood from a heparinized access line, or indwelling heparin lock;
- Discontinue fluids drip, if required;
- Use a two-syringe technique - discard the first 5 cc draw. Obtain a 3 cc sample (1 cc for P214) with the second syringe for testing.

#### Extracorporeal blood line port

- Flush the extracorporeal blood access line by withdrawing and discarding 5 cc of blood.
- Draw 3 cc sample (1 cc for P214) with a second syringe for testing.

#### Venipuncture

- Obtain a 3 cc sample (1 cc for P214) with a syringe.

*NOTE: Do not collect fresh whole blood samples using glass blood collection tubes.*

#### TEST PROCEDURE

##### FTCA510, FTK-ACT, HRFTCA510, HRFTK-ACT

- From the collection syringe, dispense exactly 2.0 cc of blood into the ACT test tube. At the same time, depress the **START** key of the appropriate test well.
- Immediately close the flip-top and agitate the test tube vigorously from end-to-end ten times.
- Insert the ACT test tube into the appropriate test well. Quickly rotate the tube clockwise. See appropriate instrument operator's manual for additional information.
- Repeat steps 1 through 3 for a duplicate ACT test tube if desired.

##### P214

- Before drawing blood, hold the test tube vertically and tap the tube on a surface to ensure glass beads are at the bottom.
- From the collection syringe, dispense exactly 0.4 cc of blood into the P214 test tube. At the same time, depress the **START** key of the appropriate test well.
- Immediately close the flip-top and flick the bottom of the tube 5-7 times with thumb and forefinger.
- Insert the test tube into the appropriate test well. Quickly rotate the tube clockwise. See appropriate instrument operator's manual for additional information.
- Repeat steps 1 through 3 for a duplicate P214 test tube if desired.
- At the sound of the buzzer, record the test result.

#### PRODUCT USE WARNINGS

*NOTE: Observe universal precautions at all times.*

- The transfer of the blood to the test tubes from a syringe should be performed using a *needleless* system. The flip-top stopper should be opened completely to dispense the blood specimen and closed completely before shaking.
- Always use a two-hand technique to transfer blood. One hand securely holds the tube while the second hand dispenses the blood specimen.
- The FTCA510, HRFTCA510, FTK-ACT and HRFTK-ACT test tubes are made of glass, which can easily be broken or cracked if mishandled. Do not drop or toss.
- The P214 test tube is a plastic tube, which can be pierced by a misplaced or a missed needle.
- Though not recommended, if a needle transfer procedure must be used, the needle must be placed in the tube parallel to the sides and centered within the tube during transfer.
- The flip top stopper should never be pierced by a needle or sharp object, due to the danger of slipping off the flip-top and piercing a finger.
- The P214 test tubes contain glass particles. Do not handle, aerosol, or ingest.
- All used test tubes containing human derived blood should be discarded in approved biohazard containers.

#### PERFORMANCE CHARACTERISTICS

**Normal Ranges**  
HEMOCHRON ACT tests were evaluated on normal, healthy volunteer donors. The expected range response is:  
**Normal Donors**

	Range:		Normal Range:
FTCA510	105-167 seconds	FTCA510	86-147 seconds
FTK-ACT	91-151 seconds	FTK-ACT	89-153 seconds
P214	110-182 seconds		

##### Precision

Reproducibility was determined using HEMOCHRON control material for ACT test tubes.

FTCA510	n	Mean (sec)	CV (%)
Normal	15	103	8.9
Ab. Level II	17	457	5.3

##### FTKACT


Normal	23	116	5.3
Ab. Level II	15	469	6.2

##### P214

Normal	24	115	6.1
Abnormal	25	228	5.2

#### Heparin Sensitivity:

Freshly obtained donor blood was heparinized *in vitro*. The graph compares heparin dose response of the ACT tests in absence of heparin.

**Figures A and B: Sensitivity of Celite and Kaolin ACT tests to Aprotinin.** Results of an *in vitro* analysis of triplicate determinations using blood samples from 2 normal donors (N=6 for each bar). Statistically significant differences between Celite ACT and kaolin ACT values are denoted as \* (p<0.05) or + (p<0.01).

*NOTE: The steady state plasma concentration of aprotinin after a full Hammersmith dose is 270 KIU/ml. As the level of aprotinin increased (600 KIU/ml), the kaolin ACT may be affected in highly heparinized blood samples.*

Fig. A. Baseline (0 heparin)



Fig. B. 3.0 units/ml heparin

#### LIMITATIONS OF THE PROCEDURE

The ACT is affected by poor technique, including both sample collection and test procedure. Proper specimen/reagent mixing is required for precise and accurate testing. The following may affect results or be misleading in test interpretation:

- Hemodilution, cardiologic solutions, hypothermia, platelet dysfunction, hypofibrinogenia, other coagulopathies, and certain medications.
- Test kits that have been improperly stored, affected by heat, or expired.

Test results which do not agree with expected values should be verified and thereafter evaluated by alternative diagnostic means.

#### QUALITY CONTROL

Routine quality control (QC) testing and tracking should be a part of a comprehensive quality assurance program. HEMOCHRON Whole Blood Coagulation System Quality Control products are available to make routine QC convenient and affordable.

#### Daily Instrument QC

At a minimum, all HEMOCHRON instruments should be quality controlled at two levels of performance, including both the normal and abnormal ranges, once every 8 hours of operation.

To assist in accomplishing daily QC, Electronic System Verification Tubes are available to provide multiple level (normal and abnormal) quality control checks on the instrument. Electronic System Verification should be performed once every 8 hours during which the instrument is utilized. This will ensure proper instrument operation.

#### QC of HEMOCHRON Test Tubes

Each box of HEMOCHRON test tubes should be validated at least once, prior to use. This can be accomplished by using the appropriate HEMOCHRON Liquid Quality Control. Acceptable performance ranges for the test tubes are included in the HEMOCHRON Quality Control Product Kit. After each individual box of test tubes has been verified, the "Performance Verified" label provided should be completed and placed on the box. This box is now "IN CONTROL" and will not require further liquid control unless a shift in clinical results is experienced.

*NOTE: If multiple boxes are received within the same shipment, it is recommended to validate each box upon opening, prior to use.*

#### REFERENCES

- Doty DB, Knot HW, Hoyt JL, Koepke JA: Heparin dose for accurate anticoagulation in cardiac surgery. J Cardiovasc Surg 1979; 20:597-604.
- Hattershy PG: Activated coagulation time of whole blood. JAMA 1966; 136-146.
- Ogilby JD, Koepelman HA, Klein LW, Agarwal JB: Adequate heparinization during PTCA: Assessment using activated clotting time. JACC 1988; 11:237A.
- Bull BS, Korpman RA, Hase WM, Briggs BD: Heparin therapy during extracorporeal circulation. I. Problems inherent in existing heparin protocols. J Thorac Cardiovasc Surg 1975; 69-674-684.
- Wang JS, Lin Cy, Hung WT, Thisted RA, Karp RB: In vitro effects of aprotinin on activated clotting time measured with different activators. J Thorac Cardiovasc Surg 1992; 104-113:5-1140.
- Sanders PW, Curtis JJ: Management of anticoagulation for hemodialysis. In: Sissenon AR, Fine RN (editors): Dialysis Therapy. The CV. Mosby Company, St. Louis, MO, 1986; pp:39-41.

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##### RÉFÉRENCES

- Dory DB, Knot HW, Hoyt JL, Koepke JA: Heparin dose for accurate anticoagulation in cardiac surgery. J Cardiovasc Surg 1979; 20:597-604.
- Hattersley PG: Activated coagulation time of whole blood. JAMA 1966; 136-146.
- Ogily JJ, Kopolman HA, Klein LW, Agarwal JB: Adequate heparinization during PTCA. Assessment using activated clotting time. JACC 1988; 11:237A.
- Bull BS, Korpman RA, Huse WM, Briggs BD: Heparin therapy during extracorporeal circulation: I. Problems inherent in existing heparin protocols. J Thorac Cardiovasc Surg 1975; 69:674-684.
- Wang JS, Lin Cy, Hung WT, Thisted RA, Karp RB: In vitro effects of aprotinin on activated clotting time measured with different activators. J Thorac Cardiovasc Surg 1992; 104: 1135-1140.
- Sanders PW, Curtis JJ: Management of anticoagulation for hemodialysis. In: Nissensoen AR, Fine RN (editors). Dialysis Therapy. The CV. Mosby Company, St. Louis, MO, 1986; pp39-41.

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