REF
R-HRT 480P
R-HRT480PK

Hemochron®

Whole Blood Coagulation System

Heparin Response Test (HRT) Package Insert English

INTENDED USE

The Hemochron® HRT is designed for *in vitro* diagnostic heparin dose response evaluation. The Hemochron HRT480 is intended for use in invasive cardiology procedures requiring moderate to high levels of systemic heparinization. The kaolin-activated HRT assays may also be used for patients receiving aprotinin. The HRT480 should be used with established Hemochron procedural guidelines with Hemochron models 401, 801, 8000 and *Response*.

For in vitro Diagnostic Use / For Professional Use / Rx Only

SUMMARY AND EXPLANATION
The anticoagulant properties of heparin are utilized during many invasive procedures to counterbalance the thrombotic tendency which results when blood is exposed to foreign surfaces. The necessary level of anticoagulation varies with the procedure as well as its duration. Frequently, the Activated Clotting Time (ACT) is performed at the bedside to monitor the anticoagulant status before, during and after invasive procedures to limit bleeding and thrombotic complications. ACT target times have been established in each setting^{1,2}, providing an indicator of adequate anticoagulation while avoiding over-heparinization. Higher heparin concentrations during bypass procedures may correlate with greater postoperative blood loss³, demonstrating the value of precise heparin regulation within the target range.

Yet even when closely monitored, effective heparin therapy is hampered by unpredictability that results from the inherent variation of different heparin preparations⁴ as well as a wide variance in patient response. Bull, et al⁵ demonstrated that the amount of heparin required to produce a target level of anticoagulation during bypass surgery varies threefold from patient to patient. Clearly, the management of heparin therapy based on empirical dosing protocols (e.g., dose/kg body weight) will be grossly inadequate to a significant number of patients⁶.

Heparin Response Test
The HRT may be used in conjunction with the injectable heparin preparation supplied by the individual hospital pharmacy. The HRT predicts appropriate heparin dosages by providing an *in vitro* evaluation of the patient's sensitivity to the heparin in the HRT tube. It determines the dose of heparin required to bring even the resistant or sensitive patient to the appropriate level of anticoagulation. Studies have shown that using the RxDx system gives optimal results and improved clinical outcomes, leading to reduction in surgical costs^{7,8}. Data show that running the HRT in this manner allows accurate prediction of the *in vivo* response to heparin on an *in vitro* assay^{9,10}.

The HRT test is available with either Celite® (diatomaceous earth) or kaolin activators. The Celite-activated HRT test is not intended for use with the protease inhibitor aprotinin (Trasylol®, Bayer Corporation), which may be administered to reduce post-operative bleeding, especially during cardiopulmonary bypass surgery, and can prolong the Celite-activated HRT. Kaolin is unaffected by moderate doses of aprotinin. The kaolin-activated HRT can be used for heparin dose response evaluation in patients whether or not they are receiving aprotinin.

The HRT is performed using Hemochron Whole Blood Coagulation instrument models 401, 801, 8000 or *Response*. The HRT contains the reagents required to perform an ACT test plus an appropriate amount of heparin from porcine intestinal mucosa. The price of the HRT is based on the amount of heparin used. Studies by Bull, et al in 1975⁵, Patient whole blood is used to perform a baseline ACT and an HRT which contains 3 units/ml of heparin. A heparin dose curve is constructed for the individual patient. Using these results and patient data, the patient's individual heparin requirement in units/ml is determined. The bolus dose is calculated based on the patient's blood volume and gender¹¹. This can be calculated automatically by the Hemochron *Response* (software version 2.0 or higher) or 8000, or manually using the worksheet provided.

REAENTS

The Hemochron HRT test tubes are intended for use with fresh whole blood. The test tubes contain a lyophilized preparation of USP characterized heparin with added activator, stabilizers and buffers.

- Each R-HRT480P tube contains the equivalent of 6 USP units of porcine mucosal heparin (3 units/ml in blood), diatomaceous earth, stabilizers and buffers.
- Each R-HRT480PK tube contains the equivalent of 6 USP units of porcine mucosal heparin (3 units/ml in blood), kaolin, stabilizers and buffers.

PRODUCT USE WARNING

NOTE: Observe universal precautions at all times.

- The blood specimen should be transferred using an appropriate transfer needle to pierce the stopper.
- Always use a two-hand technique to transfer blood. One hand securely holds the tube while the second hand dispenses the blood specimen from the syringe.
- HRT test tubes are made of glass. They can be broken or cracked if mishandled. Do not drop or toss tubes.
- The HRT test tubes contain a particulate activator (Celite or kaolin) as well as a material of biological origin (heparin). Do not handle the material, aerosol or ingest.
- All used test tubes containing human derived blood should be discarded in approved biohazard containers.
- Do not remove stoppers to deliver a blood sample to the tube.

STORAGE AND STABILITY

The Hemochron HRT test tubes should be stored at controlled room temperature (15° - 30°C). The HRT tubes are stable until the marked expiration date. Avoid prolonged exposure of the HRT kits to temperatures exceeding 30°C.

SPECIMEN COLLECTION AND HANDLING

NOTE: For open chest procedures, specimens for HRT performance should be collected post-sternotomy¹². For all other applications, specimens for HRT performance should be collected from the venous line, if possible. In all cases, to maximize accuracy of heparin prediction, the HRT should be performed within thirty (30) minutes prior to heparin infusion. If the 30 minutes is exceeded, the HRT test should be repeated.

Materials Provided

- HRT test tube (Celite or kaolin-activated as appropriate).
- HRT Worksheet (for manual calculations if required).

Materials Required

- Hemochron *Response*, 8000, 801 or two 401s.
- Non-siliconized, non-heparinized 5 mL syringes for blood specimen collection.
- #20 gauge or larger straight needle for transfer of the blood specimen into the test tube.
- ACT test tube: HRTCA510 (black top) or HRFTK-ACT (gold top), as appropriate.

CAUTION: Every precaution should be taken to use proper technique with syringes to avoid accidental needlesticks.

For blood collection, adhere to the appropriate technique (A or B):

A. Involving venous or arterial blood line. (Do not obtain blood from a heparinized access line or indwelling heparin lock.)

- Discontinue fluids drip, if required.
- Using a two syringe technique, draw two 5 mL syringes and discard the first syringe.
- Extracorporeal blood port.
 - Flush the extracorporeal blood access line by withdrawing and discarding 5 mL of blood.
 - Draw 5 mL whole blood into a 5 mL syringe.

QUALITY CONTROL (QC)

Routine quality control testing and tracking should be a part of a comprehensive quality assurance program. Hemochron Whole Blood Coagulation System Quality Control products are available to make routine QC convenient and affordable.

Daily QC of the Instrument
Hemochron instruments should be quality controlled at two levels of control every 8 hours of operation.

To assist in accomplishing daily QC, Electronic System Verification Tubes are available and can provide multiple level (normal and abnormal) quality control checks on the instrument. A Temperature Verification Tube is also provided for use every 6 months to verify that the instrument is maintaining the proper temperature 37°C ± 1°C. Any errors found in the temperature controls system will also be displayed on the LCD screen by the instrument.

QC of Hemochron Test Tubes
Each box of Hemochron test tubes should be validated at least once prior to use. This can be accomplished by using the appropriate Hemochron Liquid Quality Control. Acceptable performance ranges for the test tubes are included with each Hemochron Quality Control Product Kit. After each individual box of test tubes has been verified, the "Performance Verified" table provided on the side panel of each test tube box should be completed. This box is now "IN CONTROL" and will not require further liquid quality control unless a shift in clinical results is experienced.

NOTE: If multiple boxes are received within the same shipment, it is recommended to validate each box prior to use.

TEST PROCEDURE

If using the Hemochron *Response* (software version 2.0 or higher) or 8000, the system software will quantify the heparin dose based on the HRT test results and information supplied by the user, along with baseline ACT results. If using the Hemochron *Response* (software versions below 2.0), 401 or 801, use the HRT calculation worksheet provided with the test kit to calculate heparin dose.

If using Celite-activated HRT tests, an HRFTKA510 (black top) ACT should be used. If using kaolin-activated HRT, an HRFTK-ACT (gold top) should be used.

Recommended protocol

- Within 30 minutes prior to heparinization of the patient, simultaneously perform a baseline ACT (see ACT package insert for details) and an HRT for heparin bolus dose calculation. From the collection syringe, dispense exactly 2.0 mL of blood into the HRT tube. At the same time, depress the START button of the appropriate Hemochron test well.
- Immediately agitate the test tube vigorously from end to end ten times.
- Insert the HRT test tube into the appropriate Hemochron test well. Quickly rotate the tube clockwise. See the appropriate instrument operator manual for additional information.
- From the same collection syringe, dispense exactly 2.0 mL of blood into the ACT tube. At the same time, depress the START button of the appropriate Hemochron test well.
- Immediately agitate the test tube vigorously from end to end ten times.
- Insert the ACT test tube into the appropriate Hemochron test well. Quickly rotate the tube clockwise. See the appropriate instrument operator manual for additional information.
- At the indicator tone, record the test result.
- Determine the heparin bolus dose. See Dose Calculations below.
- Perform an ACT five (5) minutes post heparin bolus to verify that target level anticoagulation has been attained. Should a subsequent dose be required, refer to Subsequent Dosing information below.
- Monitor the ACT throughout the case, infusing a subsequent dose of heparin when necessary.

Dose Calculations

If using the Hemochron *Response* (software version 2.0 or higher), refer to the Hemochron *Response* RxDx Operator Manual for detailed instructions on the use of this automated calculation system. If using the Hemochron 8000, follow instructions in this instrument's operator manual. If using the Hemochron *Response* with software versions below 2.0, Hemochron 401 or 801, follow instructions on the HRT worksheet provided in this kit to calculate the heparin dose.

Subsequent Dosing
If target level anticoagulation is not attained or ACT results fall below target level during the procedure, a subsequent heparin dose may be determined by the Hemochron *Response* (software version 2.0 or higher) or 8000, or by applying the following formula:

Subsequent dose =

Target ACT
−
Current ACT
∗
(Post bolus ACT
−
Baseline ACT)

/

Initial Heparin Dose

*Current ACT = Post heparin ACT of pre-bypass patient / status ACT of patient on bypass

An ACT should be performed 5 minutes following infusion of the subsequent dose to verify that target level anticoagulation has been attained. As with any other *in vitro* test procedure, the HRT results should never supersede the clinical judgment of the attending physician.

LIMITATIONS OF THE PROCEDURE

The HRT is affected by poor technique in the test procedure. Proper specimen/reagent mixing is required for precise and accurate testing. The following conditions may affect test results or be misleading in test interpretation:

- Unexpected prior anticoagulation with high levels of heparin. Small quantities of heparin in the preoperative patient (baseline ACT < 220 seconds) will not affect the accuracy of the predicted dose required to reach the target time. The HRT test may be used if the baseline ACT is less than 220 seconds.
- The level of accuracy of the calculated heparin dose by any dose response method is affected by the accuracy of the blood volume calculation. Although blood volume calculations based on gender, height and weight are generally accurate, an extreme variance in expected blood volume may be seen in patients with unusually high or low amounts of adipose tissue, excessive intravenous fluid infusion as well as other medical and physiological conditions affecting blood volume.
- Other medications known to alter the patient's physiological status, (i.e., nitroglycerin, protease inhibitors) may affect the HRT results or heparin requirements. As with any other *in vitro* test procedure, the HRT test results should never supersede the clinical judgment of the attending physician.
- Only kaolin-activated R-HRT480PK tubes should be used for patients receiving aprotinin. Studies have not been performed utilizing other antifibrinolytic agents.
- In vitro* analyses have indicated an elevation of the HRT predicted dose at aprotinin concentrations exceeding 500 KIU/mL blood when using the kaolin-activated system.¹³
- Test results which do not agree with expected values should be verified and thereafter evaluated by alternative diagnostic means.

PERFORMANCE CHARACTERISTICS

The HRT result is an indication of the heparin sensitivity of the individual patient. Lower HRT values indicate a more heparin resistant patient requiring a higher heparin bolus dose to reach target anticoagulation levels. Published clinical reports have shown that HRT recommended doses can vary from 50% to 200% of standard empirical dosing.¹⁴

	Level 1				Level 2			
R-HRT480P	Day 1	Day 2	Day 3	Combined	Day 1	Day 2	Day 3	Combined
N	3	3	5	11	3	3	5	11
Mean	137	140	136	136	337	344	343	342
Std dev	1.0	7.0	9.0	6.1	2.6	6.7	8.7	7.1
CV%	0.7	5.0	6.6	4.5	0.8	1.9	2.5	2.1

R-HRT480PK	Day 1	Day 2	Day 3	Combined
N	3	3	5	11
Mean	139	134	151	134
Std dev	0.6	2.9	1.0	3.4
CV%	0.4	2.2	0.8	2.5
	3.0	3.8	0.3	2.7

Kaolin-activated dosing results in cardiac surgery patients treated with aprotinin
Clinical data were collected at two sites following institutional approval. Dosing calculations at each site were performed manually by the perfusionist running the assays.

The efficacy of the HRT predicted doses are summarized in the table below. The percentage of patients given the dose as quantified by the kaolin HRT test is shown under the heading "Dosed as recommended"

Heparin Dosing

# Patients	Dosed as recommended	Target ACT	Pts Reaching Target	Heparin Dose (x 1000IU/L)
Site	#	(%)	#	(%)

1	77	71	(92)	400	66	(93)	Not reported
2	42	39	(93)	450	38	(97)	287 ± 10.6 12.5 - 68.0

These data reflect dosing efficacy of the kaolin system in the presence of aprotinin. Wide dose ranges are expected in a population study due to the individual patient's heparin sensitivity.

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