HEMOCHRON®

Whole Blood Coagulation Systems

HITRx Whole Blood Controls

Package Insert

INTENDED USE
The HEMOCHRON® HITRx Whole Blood Quality Controls are intended to be used to perform quality control assays using the HEMOCHRON HTT and HEMOCHRON RxDs™ kaolin and Celite® activated HRT480, PRT400 and PDAO test tubes.

SUMMARY AND EXPLANATION

Blood coagulation instruments and assays should be quality controlled prior to and during routine use. A proper quality assurance program including regular instrument maintenance, quality control assay measurements and complete procedure records.

REAGENTS
Each box of HEMOCHRON HITRx Whole Blood Quality Controls contains both Level 1 and Level 2 quality control materials. Assayed clotting times are provided with each lot of material.

- Each HITRx quality control product consists of:
 - Dried fixed bovine red blood cells, buffered sheep and horse plasma in individual test vials.
 - Diluent preparation, consisting of distilled water, sodium chloride, Tween 20, calcium chloride, anticoagulant and preservatives in individual vials for reconstitution of the whole blood control. HRT quality control diluents also contain protamine sulfate.

NOTE: Diluents are lot matched to the whole blood control vials. DO NOT use diluents from one lot of control with whole blood vials from a different lot.

Three Level 1 and Level 2 quality control products are available for use with the following HEMOCHRON tests in the HEMOCHRON Model 401, 801, 8000 and Response instruments.

- Level 1 and Level 2 quality control, cat. # RQCHRT: for use with HEMOCHRON kaolin and Celite activated bovine or porcine HITRx tubes.
- Level 1 and Level 2 quality control, cat. # RQCPRT: for use with HEMOCHRON kaolin and Celite activated PRT200, PRT400 and PDAO tubes.
- Level 1 and Level 2 quality control, cat. # RQCHTT: for use with HEMOCHRON HTT tubes.

CAUTION: Animal blood products are included in the HITRx quality control preparation. No material in any preparation is derived from human blood. All blood products, whether of human origin or animal origin, should be handled with care.

HANDLING AND STORAGE
When refrigerated (2-8°C), the vials are stable until the marked expiration date. The quality control product should never be exposed to temperatures in excess of 37°C. Reconstituted vials should be used immediately. HITRx Quality Controls may also be stored at room temperature for up to 4 weeks. (The marked expiration date must not be exceeded.) If the product is to be stored at room temperature, cross out the stamped expiration date on the right side of the box, and fill out the redating table printed on the front panel of the box with the new expiration date. Any change should be properly initialed and dated.

PROCEDURE
Materials Provided
Quality control product kits for use with the HEMOCHRON HRT480, PRT200, PRT400, PDAO and HTT tubes contain the following:
• 10 vials of dried Level 1 whole blood control (2.0 cc).
• 10 vials of Level 1 diluent (4.0 cc).
• 10 vials of dried Level 2 whole blood control (2.0 cc).
• 10 vials of Level 2 diluent (4.0 cc).
• A quality control tracking chart.

Materials Required But Not Provided
• HEMOCHRON Model 401, 801, 8000 or Response instrument.
• HEMOCHRON HRT480, PRT200, PRT400, PDAO or HTT tubes.
• 3 cc plastic syringe with a 20 or 21 gauge needle for control vial and HITT tube rehydration.

Reagent Preparation
A. Quality Control Product
Remove the appropriate quality control vials from the refrigerator and allow them to come to room temperature (15-30°C) prior to testing. This could require up to 60 minutes.

B. HEMOCHRON test tubes
Remove the test tubes from the refrigerator and allow them to come to room temperature prior to testing. This could require up to 60 minutes.

CAUTION: All used tubes and control vials should be considered as potentially infectious, handled with care, and disposed of in accordance with standard medical waste disposal policy.

IF OUT-OF-RANGE RESULTS OR FAULT MESSAGES PERSIST, CONTACT ITC TECHNICAL SERVICES.

TROUBLESHOOTING
Please refer to the appropriate instrument Operator's Manual "Troubleshooting" section for instrument software fault messages, definitions and actions.

PROBLEM	CAUSE	ACTION
Control value is above or below published range	Reconstituted control was not thoroughly mixed.	Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.
	Time period between control material mixing and addition to test tube is too long	Repeat test. After thorough mixing of the control material and diluent, IMMEDIATELY withdraw the rehydrated whole blood control and transfer to the appropriate HEMOCHRON tube.
Fault message was obtained (see appropriate instrument Operator's manual)	Reconstituted control was not thoroughly mixed.	Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

L1	L2	\bar{x}	r
EN Level 1	Level 2	Mean	Range
DE Level 1	Level 2	Durchschnitt	Bereich
FR Niveau 1	Niveau 2	Moyenne	Plage
IT Livello 1	Livello 2	Media	Intervallo
ES Nivel 1	Nivel 2	Media	Intervalo
PT Nivel 1	Nível 2	Media	Limite
NK Niveau 1	Niveau 2	Bereik	Bereik
DK Niveau 1	Niveau 2	Gennemsnit	Område

REFERENCES
1. Statland, B and Westgard, J: Quality Control: Theory and Practice, "In Clinical Diagnosis and Management by Laboratory Methods", ed, Henry, J (Philadelphia WB Saunders Co., 1984) pp.74-93.
2. Colton, C: Statistics in Medicine (Boston: Little, Brown, and Co., 1974) pp 289-313.

HEMOCHRON® is a registered trademark of ITC. Celite® is a registered trademark of Celite Corporation.

8 Olsen Avenue • Edison, NJ 08820 USA
Tel: 732.548.7700 • Fax: 732.248.1928
www.itctmed.com

a subsidiary of Thoratec Corporation

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PROBLEM CAUSE ACTION
Control value is above or below published range

Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

Repeat test. After thorough mixing of the control material and diluent, IMMEDIATELY withdraw the rehydrated whole blood control and transfer to the appropriate HEMOCHRON tube.

Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

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