REF RQCHRT

REF RQCPRT

Hemochron®

Whole Blood Coagulation System

HiTRx Whole Blood Controls

Package Insert

INTENDED USE

The HEMOCHRON® HiTRx Whole Blood Quality Controls are intended to be used to perform quality control assays using the HEMOCHRON HiTT and HEMOCHRON RxDx® kaolin and Celite® activated HRT480, PRT200, PRT400 and PDAO test tubes.

For *in vitro* Diagnostic Use

SUMMARY AND EXPLANATION

Blood coagulation instruments and assays should be quality controlled prior to and during routine use. A proper quality assurance program includes regular instrument maintenance, quality control assay measurements and complete performance records.

REAGENTS

Each box of HEMOCHRON HiTRx Whole Blood Quality Controls contains both Level 1 and Level 2 quality control materials. Assayed clotting times are provided with each lot of material.

Each HiTRx quality control product consists of:

- Dried fixed bovine red blood cells, buffered sheep and horse plasma in individual test vials.
- Diluent preparation, consisting of distilled water, sodium chloride, Tween 20, calcium chloride, anticoagulant and preservatives in individual vials for reconstitution of the whole blood control. HRT quality control diluents also contain protamine sulfate.

NOTE: Diluents are lot matched to the whole blood control vials. DO NOT use diluents from one lot of control with whole blood vials from a different lot.

Three Level 1 and Level 2 quality control products are available for use with the following HEMOCHRON tests in the HEMOCHRON Model 401, 801, 8000 and Response instruments.

- Level 1 and Level 2 quality control, cat. # RQCHRT: for use with HEMOCHRON kaolin and Celite activated bovine or porcine HRT480 tubes.
- Level 1 and Level 2 quality control, cat. # RQCPRT: for use with HEMOCHRON kaolin and Celite activated PRT200, PRT400 and PDAO tubes.
- Level 1 and Level 2 quality control, cat. # RQCHITT: for use with HEMOCHRON HiTT tubes.

CAUTION: Animal blood products are included in the HiTRx quality control preparation. No material in any preparation is derived from human blood. All blood products, whether of human origin or animal origin, should be handled with care.

HANDLING AND STORAGE

When refrigerated (2-8 °C), the vials are stable until the marked expiration date. The quality control product should never be exposed to temperatures in excess of 37 °C. Reconstituted vials should be used immediately. HiTRx Quality Controls may also be stored at room temperature for up to 4 weeks. (The marked expiration date must not be exceeded.) If the product is to be stored at room temperature, cross out the stamped expiration date on the right side of the box, and fill out the redating table printed on the front panel of the box with the new expiration date. Any change should be properly initialed and dated.

PROCEDURE

Materials Provided

Quality control product kits for use with the HEMOCHRON HRT480, PRT200, PRT400, PDAO and HiTT tubes contain the following:

- 10 vials of dried Level 1 whole blood control (2.0 cc).
- 10 vials of Level 1 diluent (4.0 cc).
- 10 vials of dried Level 2 whole blood control (2.0 cc).
- 10 vials of Level 2 diluent (4.0 cc).
- A quality control tracking chart.

Materials Required But Not Provided

- HEMOCHRON Model 401, 801, 8000 or Response instrument.
- HEMOCHRON HRT480, PRT200, PRT400, PDAO or HiTT tubes.
- 3 cc plastic syringe with a 20 or 21 gauge needle for control vial and HiTT tube rehydration.

Reagent Preparation

A. Quality Control Product
Remove the appropriate quality control vials from the refrigerator and allow them to come to room temperature (15-30°C) prior to testing. This could require up to 60 minutes.

B. HEMOCHRON test tubes

Remove the test tubes from the refrigerator and allow them to come to room temperature prior to testing. This could require up to 60 minutes.

CAUTION: All used tubes and control vials should be considered as potentially infectious, handled with care, and disposed of in accordance with standard medical waste disposal policy.

NOTE: Reconstitution and mixing of the whole blood material should be accomplished quickly and without delay in any step. Once the dried material has been reconstituted, the sample should be used immediately, as clotting will occur.

	L1	L2	\bar{x}	r
EN	Level 1	Level 2	Mean	Range
DE	Level 1	Level 2	Durchschnitt	Bereich
FR	Niveau 1	Niveau 2	Moyenne	Plage
IT	Livello 1	Livello 2	Media	Intervallo
ES	Nivel 1	Nivel 2	Media	Intervalo
PT	Nivel 1	Nivel 2	Media	Limite

HEMOCHRON® is a registered trademark of Accriva. Celite® is a registered trademark of Celite Corporation.

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Test procedure for HEMOCHRON HRT, PRT and PDAO

- Peel back the flap of the crimp seal on the diluent vial ensuring that the Level on the diluent vial matches the Level on the Whole Blood QC vial, but do not remove the entire crimp seal (use this technique when opening each vial). Using a 3 cc syringe with a needle, withdraw exactly 3.0 cc from the diluent vial. Transfer the diluent into the whole blood vial by direct puncture of the stopper.
- Remove the syringe and needle from the whole blood vial. Do not discard syringe and needle.
- Using moderate end-to-end inversion, mix the control vial for EXACTLY 15 seconds.
- Using the same syringe, withdraw 2.0 cc of the rehydrated whole blood control.
- Immediately transfer the 2.0 cc into the appropriate HEMOCHRON tube. As the control material first enters the tube, depress the START button on the appropriate HEMOCHRON test well.
- Immediately agitate the tube vigorously from end to end 10 times.
- Insert the tube into the appropriate HEMOCHRON test well. Quickly rotate the tube clockwise. See the appropriate instrument operator manual for additional information.
- Upon test completion, record results.

Test procedure for HEMOCHRON HiTT

Preparation of the HiTT tube

- Peel back the flap of the crimp seal on the distilled water vial found in the HiTT test kit, but do not remove the entire crimp seal (use this technique when opening each vial). Using a 3 cc syringe with a needle, withdraw exactly 0.5 cc from the distilled water vial. Transfer the water into the turquoise-top HiTT tube by direct puncture of the stopper.
- Shake the HiTT tube gently from end-to-end until the lyophilized activator is thoroughly dissolved.

NOTE: For optimal performance, HiTT tubes should be incubated within thirty (30) minutes after rehydration.

- Insert the HiTT test tube into the appropriate HEMOCHRON test well. To begin the timer for the 180 second HiTT test incubation period, depress the START button or the PRE-WARM/SELECT button for the appropriate test well. For the Response instrument, press the MENU key of the appropriate test well. Press the #3 key for PreWarm Well. Then press the #3 key to select 180 – All TT; this will begin the prewarm function.
- After 180 seconds have elapsed on the HEMOCHRON instrument, remove the HiTT tube from the instrument.

NOTE: It is important to remove the HiTT test tube from the HEMOCHRON instrument to avoid over-incubation.

Preparation of the Quality Control

- Peel back the flap of the crimp seal on the diluent vial, but do not remove the entire crimp seal. Using a 3 cc syringe with a needle, withdraw exactly 3.0 cc from the diluent vial this can be the same syringe used to reconstitute the HiTT tube). Transfer the diluent into the whole blood vial by direct puncture of the stopper.
- Remove the syringe and needle from the whole blood vial. Do not discard syringe and needle.
- Using moderate end-to-end inversion, mix the control vial for EXACTLY 15 seconds.
- Using the same syringe, withdraw 1.5 cc of the rehydrated whole blood control.
- Immediately transfer the 1.5 cc into the HiTT tube. As the control material first enters the tube, depress the START button on the appropriate HEMOCHRON test well.
- Immediately agitate the tube vigorously from end to end 10 times.
- Insert the tube into the appropriate HEMOCHRON test well. Quickly rotate the tube clockwise. See the appropriate instrument operator manual for additional information.
- Upon test completion, record results.

PERFORMANCE CHARACTERISTICS

Expected Values

The acceptable range for each Quality Control Product is published on the back of this package insert. These values have been established at the manufacturer's facility and are lot specific. The range represents the acceptable performance range for all HEMOCHRON tube instruments. The user should be able to obtain a value within the published range included for each product. A proper quality assurance program includes regular control measurements and maintenance of complete performance records.

Each institution may also establish its own expected range of response based on the mean and standard deviation of at least 20 repeated tests. The institution's range should be narrower and fall within the published range listed on this package insert. Studies show that intra-laboratory variation in test results should produce a coefficient of variation of approximately 10% or less for coagulation control tests.

LIMITATIONS

Technique variation can significantly alter performance characteristics. Care must be taken to establish and adhere to the recommended test procedure found in this package insert and the associated test tube package insert. Variations in ambient temperature should be minimized for consistency of test performance. In cases where results fall outside of the assigned range, patient results should be considered suspect and performance of quality control tests should be repeated. Each institution should establish a data tracking record to identify instruments, tubes, reagents or operators whose performance is outside of the institution's expected range.

In cases where results are outside of the expected range, the cause will likely fall into one of the following possible categories:

- Operator technique
- Quality Control Product/Test Tube
- Instrumentation
- Environmental factors

IF OUT-OF-RANGE RESULTS OR FAULT MESSAGES PERSIST, CONTACT ITC TECHNICAL SERVICES.

TROUBLESHOOTING

Please refer to the appropriate instrument Operator's Manual "Troubleshooting" section for instrument software fault messages, definitions and actions.

PROBLEM	CAUSE	ACTION
Control value is above or below published range	Reconstituted control was not thoroughly mixed.	Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.
	Time period between control material mixing addition to test tube is too long	Repeat test. After thorough mixing of the control material and diluent, IMMEDIATELY withdraw the rehydrated whole blood control and transfer to the appropriate HEMOCHRON tube.
Fault message was obtained (see appropriate instrument Operator's manual)	Reconstituted control was not thoroughly mixed.	Repeat test using moderate end-to-end inversion, and thoroughly mix the control vial for EXACTLY 15 seconds.

REFERENCES

- Statland, B and Westgard, J: Quality Control: Theory and Practice, "In Clinical Diagnosis and Management by Laboratory Methods", ed, Henry, J (Philadelphia WB Saunders Co., 1984) pp.74-93.
- Colton, C: Statistics in Medicine (Boston: Little, Brown, and Co., 1974) pp 289-313.

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