Surgicutt® Bleeding Time Device

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Summary and Explanation
When vessels are injured, platelets adhere to the exposed subendothelial lining of the vessel and exposed collagen. Following this adhesion phase, platelets begin to aggregate at the wound site forming the primary platelet plug. The purpose of the bleeding time test is to provide a measure of such platelet functions in small vessel hemostasis¹.

Principle of the Procedure

The test is performed by making a minor standardized incision on the forearm, then wicking the blood away from the wound site using filter paper. Care must be taken to insure that the filter paper does not touch the incision during the wicking process so as not to disturb the formation of the platelet plug.

A prolonged bleeding time, based on a laboratory's normal range, does not in itself diagnose underlying platelet disorders, either qualitatively or quantitatively. It does, however, indicate the need for more quantifiable testing, such as platelet count, platelet aggregation studies, etc.

Surgicutt® is a sterile, standardized, easy-to-use, disposable device that makes a uniform, surgical incision, 5 mm long by 1 mm deep.

A spring-activated surgical steel blade is housed in a plastic unit from which it protracts and retracts automatically. Swift blade action and a standardized incision depth well above the pain fibers virtually eliminate pain sensation¹.

Surgicutt has a completely flat blade slot surface, enabling it to be placed flush against the patient's forearm without exerting undue pressure or indenting the skin, thereby assuring incisions of uniform depth and length.

Surgicutt is a single-incision device.

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- Fully-automated incision-making device
- Disposable
- Sterile
- Easy to use
- Standardized surgical-type incision
- Retractable blade
- Reattachable trigger safety clip
- Latex free

	Specifications	
	Depth	Width
Surgicutt	1.0mm	5.0mm
<i>Surgicutt devices are for single-use.</i>		

Test Procedures

Caution: There is occasional scarring when the bleeding time test is administered¹. Patients, parents or guardians should be advised of this possibility, especially in cases of a positive clinical history of keloid scar formation. Potential scarring may be reduced by approximating the skin edges with a non-allergenic wound closure strip for 24 hours¹. If there is ooze from the incision as may be encountered in severe primary hemostatic disorders, then a pressure-type dressing should be used in conjunction with the wound closure strip².

Materials Provided

- Each blister pack contains one Surgicutt device. Product is sterile unless package has been damaged or broken. If blister pack has been opened, do not use.

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DEUTSCH / GERMAN

Materials Required (not provided)

- Surgicutt Bleeding Time Blotting Paper
- Antiseptic swab
- Sphygmomanometer
- Stopwatch (with sweep second hand)
- Wound closure strip

Selecting the Site and Direction of Incision

Bleeding time test incisions are best performed over the lateral aspect, volar surface of forearm, approximately 5 cm below the antecubital crease¹. Median incisions tend to be more painful, are less reproducible, and have a higher incidence of scar formation¹. The horizontal incision with (40 mm Hg) venostasis is the most sensitive technique.

Procedure

1. Place the patient's arm on a steady support with the volar surface exposed. Lightly shave the area if body hair will interfere with the test.

2. Place the sphygmomanometer cuff on the upper arm. Inflate the cuff to 40 mm Hg.

The time between inflation of cuff and incision should be 30-60 seconds. Hold at this exact pressure for the duration of the test.

3. Cleanse the area with an antiseptic swab and allow to air dry.

4. Remove Surgicutt from the blister pack, being careful not to contaminate or touch the blade-slot surface.

5. Remove safety clip (safety clip may be replaced if the test is momentarily ,delayed, however, prolonged exposure of Surgicutt to uncontrolled environmental conditions prior to use may affect its sterility). Once the safety clip is removed, DO NOT push the trigger.

6. Hold the device securely between the thumb and the middle finger.

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7. Gently rest the Surgicutt device on the patient's forearm and apply minimal pressure so that both ends of the device are lightly touching the skin. The incision is best performed over the lateral aspect, volar surface of the forearm, parallel to and 5 cm below the antecubital crease. Avoid surface veins, scars, bruises and edematous areas. This technique is the most sensitive.

8. Gently push the trigger, starting the stopwatch simultaneously.

9. Remove Surgicutt from the patient's forearm immediately after triggering.

10. Every 30 seconds, wick the blood with blotting paper. Bring the Surgicutt blotting paper close to the incision, but DO NOT touch the paper directly to the incision, so as not to disturb the formation of a platelet plug.

11. Wick the blood every 30 seconds until blood no longer stains the paper. Stop the timer. Bleeding time is determined to the nearest 30 seconds.

12. Remove the sphygmomanometer cuff and cleanse the incision site with an ,antiseptic swab. Potential scarring can be reduced by approximating the skin edges with a non-allergenic wound closure strip2 for 24 hours¹.

! Used devices are potentially infectious. Handle with appropriate care and dispose of used devices in accordance with standard methods of biohazard control.

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Normal Population Values

Normal range studies were conducted for the Surgicutt device resulting in normal ranges of:

Surgicutt 2.0 min - 8.0 min.

It is recommended that each institution establish its own normal population range due to laboratory and population variables.

Aspirin Tolerance Test^{4,5}

Aspirin (acetylsalicylic acid) affects normal platelet function through a complex chemical mechanism. This hemostatic change may be demonstrated utilizing a sensitive bleeding time test.

The sensitivity of the Surgicutt bleeding time device to aspirin induced platelet dysfunction was studied in healthy, medication free male and female subjects through bleeding time testing before and after aspirin ingestion.

	Bleeding Time Device	Aspirin Dosage	Pop. Size	Mean	Standard Deviation
	Surgicutt pre-aspirin	---	21	4.3 min	1.2 min
	Surgicutt post-aspirin (4 hrs later)	975 mg	21	7.7 min	3.6 min

Prolongation of the Bleeding Time from Pre-aspirin to Post-aspirin Test Results

	Delta Bleeding Time	Standard Deviation of Delta
Surgicutt	3.5 minutes	2.9 minutes

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Conditions Affecting Test Results

Since the bleeding time is a screening test for both quantitative and qualitative platelet function, a prolonged bleeding time may be indicative of the following conditions:

- Thrombocytopenia
- Von Willebrand's disease and von Willebrand's variants
- Thrombasthenia
- Congenital or acquired coagulation deficiencies
- Storage Pool disease
- Bernard Soulier syndrome
- Uremia
- Liver failure