DEUTSCH / GERMAN

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Principles of the Tenderlett®

Tenderlett is a safe automated skin incision device for fingertip blood sampling. Each Tenderlett device is engineered to incise to the least invasive yet effective depth that will reach the capillaries and small vessels in the finger. They have been carefully calibrated to substantially reduce pain and improve fingertip blood sampling. Optimal blood flow is achieved and sample contamination and hemolysis are minimized. When the device is triggered, a surgical steel blade swiftly protracts at an angle of 30 degrees, produces an incision in only the uppermost vasculature of the finger, and then automatically retracts. The permanently retracting blade eliminates danger and contamination hazards. The swift blade action, produces optimal blood flow at a much safer and more comfortable wound depth. Good blood flow also minimizes sample contamination and hemolysis. The surgical blade incises and retracts rapidly resulting in a minor, superficial wound that is virtually pain-free.

Specifications

	Incision Depth	Incision Length
Tenderlett	1.75 mm	0.94 mm
Sterile, single use, latex free.		

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Special Consideration

Incision Depth

For older, normally developed, robust children, teenagers, and adults, the optimal and safe finger incision depth is 1.75 mm. Regardless of the patient’s age, great care is required for safe fingerstick blood sampling.

Incision Test Procedures

Limitations

Make sure that the fingertip skin incision site is not swollen or edematous. Do not collect blood from recent incision sites or open cuts in the finger. Do not use povidone iodine to disinfect skin incision sites due to the potential for elevated chemistry values. Cold fingers, must be pre-warmed prior to performing a fingerstick to ensure good results. It is advisable to prewarm fingers to 41 degrees centigrade when collecting specimens for blood gases. Immunosupressed patients are susceptible to infection from a fingerstick. Those patients should be under the careful observation and supervision of a physician.

Materials Provided

- Each packet contains one Tenderlett device. Each device is sterile unless package has been damaged or opened

Materials Required but not Provided

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- Blood collection container, e.g. capillary tubes
- Fresh gloves
- Sharps disposal unit
- Capillary tube sealer, if needed

Site Selection (all ages and fingers)

Correct incision site selection is vital to the success of the blood collection procedure. The middle and fourth fingers are preferred for virtually pain-free skin incisions because they are not as sensitive to pain as the index finger, and they bleed better. The index finger should be avoided if possible. Also avoid the middle portion of any finger because of increased sensitivity to pain and less robust bleeding. (See Figure 1) The incision should be made perpendicular to (across) the fingerprint approximately 3.00 mm to 5.0 mm from the nail bed, at the point where the finger begins to curve. (See Figure 2) This area of the finger is more vascularized and less sensitive to pain.

- Seat the ambulatory patient in a standard phlebotomy chair with an arm board. This position improves fingertip blood pressure and perfusion. When possible, the torso of supine patients should be raised slightly and/or the arm lowered to a position slightly below the level of the breastbone.
- Gently massage the entire length of the finger to increase the temperature and improve perfusion. (See Figure 3)

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
- Clean the incision site and surrounding area with 75% isopropanol. Thoroughly dry the site with a sterile gauze pad to prevent rapid hemolysis caused by residual alcohol.

- Remove the Tenderlett device from its plastic packet, taking care not to touch or contaminate the blade-slot surface or contoured end of the device. (See Figure 4)
- Again, gently massage the lower portion of the finger while avoiding the fingertip incision site. Firmly grasp the lower portion of the finger to restrict return circulation. Firmly position the Tenderlett device at the incision site and depress the trigger.(See Figure 5)
- After triggering, immediately remove the device from the patient’s finger. Using a sterile gauze pad, gently wipe away the first small droplet of blood. Apply gentle, continuous pressure to the finger avoiding excessive massaging as this may contaminate the sample or cause hemolysis.
- Taking care not to make direct contact between the wound and the collection container, fill to the desired specimen volume. (See Figure 6)
- Following blood collection, gently press a dry sterile gauze pad to the incision site until bleeding stops. If indicated, bandage the finger.

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

If package containing device is damaged or open, do not use.

Used devices are potentially infectious. Handle with appropriate care and dispose of used devices in accordance with standard methods of biohazard control.

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Tenderlett® is a registered trademark of Accriva Diagnostics, Inc. Patent# 5,133,730 & 5,212,879

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References

- Barr ML, Kiernan JA: The Human Nervous System. An Anatomical Viewpoint. Philadelphia: JB Lippincott Co., 1988, 39.
- Blumenfeld TA, Hertelendy WG, Ford SH: Simultaneously obtained skin puncture serum, skin puncture plasma and venous serum compared and effects of warming the skin before puncture. Clin Chem 1977;23:705-1710.
- College of American Pathologists: So you’re going to collect a blood specimen: An Introduction to Phlebotomy, 1989.
- Daly CH: The Biomechanical Characteristics of Human Skin. Thesis, University of Strathclyde, Scotland, 1966.
- Grazaitis DM, Sexson WR: Erroneously high dextrostix values caused by isopropyl alcohol. Pediatrics 1980;66:221-223.
- Henry JB: Clinical Diagnosis and Management by Laboratory Methods. Philadelphia: W.B. Saunders, 1984;3.
- Ladenson JH: Non-analytical sources of laboratory error in pH and blood gas analysis. ABSTRACTS National Bureau of Standards, June 1977.
- Landis EM: Microinjection studies of capillary blood pressure in human skin. Heart, 1930;15:209.
- Latman NS: Evaluation of electronic, digital blood glucose monitors. Biomedical Instrumentation & Technology 1991;25:43-49.
- Meites S: Comment on NCCLS guidelines for skin-puncture. Clin Chem 1986;32:231.
- Meites S, Glassco, KM: Studies on the quality of specimens obtained by skin-puncture of children II. An analysis of blood-collecting practices in a pediatric hospital. Clin Chem 1985;31:1669-1672.

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- Meites S, Levitt MJ: Skin-puncture and blood-collecting techniques for infants. Clin Chem 1979;25:183-189.
- Moritz AR, Henriques FC Jr: Studies of thermal injury, II. The relative importance of time and surface temperature in the causation of cutaneous burns. Am J Clin Pathol 1947;23:1705-1710.
- National Committee for Clinical Laboratory Standards: Procedures for the collection of blood specimens by skin punctures; approved standard. NCCLS publication H4-A2. Villanova, PA.: NCCLS ; 1986.
- O’Brien D, Ibbott FA, Rodgerson DO: Laboratory Manual of Pediatric Micro Biochemical Techniques. New York: Hoeber, 1968.
- Perry AW, Gottlieb LJ, Krizek TJ, Zachary LS: Fingerstick felons. Plastic and reconstructive surgery. Section of Plastic and Reconstructive Surgery. Chicago: Little, Brown and Company, 1988.
- Reiner CR, Meites S, Hayes JR: Optimal sites and depths for skin puncture of infants and children from anatomical measurements. Clin Chem 1990;39:547-549.
- Van Steirtegham AC , Young DS: Providone-iodine (Beta-dine) disinfectant as a source of error. Clin Chem 1977;23:1512.
- Wilkinson RH: Chemical Micromethods in Clinical Medicine, Springfield, IL: Charles C. Thomas, 1960 19-25.

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