



Incision Devices

Discover which is best for your patient's needs



Tenderfoot®

is a safe, effective device for obtaining blood samples from the heels of infants.

- The result is a high quality blood sample with significantly less pain and trauma
- The Tenderfoot produces a consistent incision that is controlled to a depth just above the nerve fibers
- Available in four different incision depths from micro-preemie to toddler



Tenderlett®

finger incision devices are designed to fit the contour of the finger.

- Available in a standard size to accommodate older children to adults
- A uniform length and depth with every incision procedure
- Safety trigger allows the blade to retract without potential harm to patient or user



Surgicutt®

bleeding time devices are used to screen for platelet function.



- All-in-one bleed time test uses the Ivy Method to screen for platelet function
- Individually packaged bleeding time blotting paper is available to wick away residual blood
- Available in a standard size to accommodate older children to adults

Softsweep® Technology

Our patented softsweep devices perform a gentle, sweeping arc incision above the level of sensitive nerve fibers.

The blade permanently retracts back into the device, eliminating any chance of an accidental stick and creates a high quality incision resulting in a sample with maximum blood flow to reduce the need for re-sticking.

Our devices are safe, easy to use, and virtually pain free with less bruising as compared with traditional lancet devices.



Illustrated above: Accriva's softsweep incision cross section (left) versus stabbing motion of traditional lancet (right).

Ordering Information

Catalog Number	Description	Depth	Length	Color	Indications*
TFM50I	Tenderfoot Micro-Preemie	0.65 mm	1.40 mm	Blue	<1000 grams
TFP50I TFP200I TFP1000I	Tenderfoot Preemie	0.85 mm	1.75 mm	White	low birth weight 1000g–2500g
TF50I TF200I TF1000I	Tenderfoot Newborn	1.00 mm	2.50 mm	Pink/Blue	birth to 6 months 2500g–9kg
TFT50I	Tenderfoot Toddler	2.00 mm	3.00 mm	Pink	6 months to 2 years >9kg
TL100I TL5000I	Tenderlett	1.75 mm	0.94 mm	Red Trigger	Adults/ older children
SU50I	Surgicutt Adult	1.00 mm	5.00 mm	Gold/White	Adult
SUP100I	Blotting Paper	9.0 cm	N/A	White	All ages

* These are suggested indications; ultimately the medical staff's clinical judgement should determine which device to use.

Tenderfoot®

Tenderlett®

Surgicutt®



To learn more, call 1-858-263-2501 or visit www.accriva.com

MINC5011 00