

Hemochron[®]
RESPONSE

Whole Blood Coagulation System

Hemochron ACTs Set The Standard

For all your hemostasis management needs



Confidently monitor coagulation at point of care

Accriva
diagnostics

Accurate Results in Any POC Setting

For more than 30 years, Hemochron® has been the standard for performing whole blood coagulation monitoring at the point-of-care.

Cardiovascular Surgery

CARDIOPULMONARY BYPASS PROCEDURES

- Heparin and Protamine System
 - Estimates the heparin doses to reach target time
 - Estimates the protamine dosage to neutralize heparin to baseline ACT result
- Kaolin based ACTs (FTK-ACT)
 - Designed for patients on aprotinin

Cardiac Catheterization Laboratory

- Reliable sheath removal in PTCA

Critical Care

DIALYSIS / ECMO


- Glass Bead Activator
 - Designed especially for low dose heparin monitoring
- Hemochron ACTs are the product of choice in most ECMO Centers

Electronic System Verification Tube

ELECTRONIC QUALITY CONTROL

- 3 Levels of EQC in one
- Reusable cost effective





CONVENIENT HANDLE
Portable bedside testing


PRINTOUT
Document and verify results

NUMERIC KEYPAD
Input patient I.D., input operator I.D.

DUAL-WELL SYSTEM
Speciality tests, duplicate tests

RS232 SERIAL PORT
Downloading

- **SPLIT SCREEN**
Easy to read
- **BAR CODE TUBES**
Positive sample identification
- **DATA MANAGEMENT**
Patient database storage
- **BATTERY OR AC OPERATED**



Ordering Information

Catalog Number	Description
HRS220	Hemochron Response Whole Blood Microcoagulation System (international)
HRS110	Hemochron Response Whole Blood Microcoagulation System
RXDX®	Analysis module for Heparin and Protamine dosing
R-HRT480P	Heparin Response Test - Porcine
R-HRT480PK	Heparin Response Test, same as R-HRT480P, but kaolin based
R-PRT200	Protamine Response Test
R-PRT200K	Protamine Response Test, same as R-PRT200, but kaolin based
R-PRT400	Protamine Response Test
R-PRT400K	Protamine Response Test, same as R-PRT400, but kaolin based
PDAO	Protamine Dose Assay
PDAOK	Protamine Dose Assay, same as PDAO, but kaolin based
HRFTCA510	Diatomaceous earth activated, non-evacuated glass test tubes
HRFTK-ACT	Kaolin activated, non-evacuated glass test tubes
P214	Glass-particulate activated, non-evacuated plastic test tubes
HE-ESV	Electronic Quality Control (EQC)

Find out More

Hemochron Response is the proven point-of-care coagulation monitoring system. Ideal for the wide variety of applications found in today's busy health care facilities.

To learn more, visit www.accriva.com

US Customers: 800.643.1640 | International Customers: 858.263.2501



Accriva
diagnostics

Accriva Diagnostics, Inc. | 6260 Sequence Drive | San Diego, CA 92121 USA
Tel: 858.263.2300 | Toll-Free: 800.631.5945 | www.accriva.com